Industrial[™] Flanktin

PBR 100L Installation Requirements





Customer is responsible for ensuring **all** the following requirements are met. Contact support@industrialplankton.com with any questions or concerns.

Consumables

pH probe solutions
☐ 7 calibration solution
☐ 10 calibration solution
☐ Storage solutions
☐ Deionized water (for rinsing)
Cleaning solutions
 ☐ 70% alcohol (For any of rubbing, ethyl, or isopropyl alcohol) ☐ Biofilm removal agents (see User Manual for alternatives) ○ Muriatic Acid (33% HCl) ○ Steris CIP 100® ☐ Bleach (4-12% sodium hypochlorite) ☐ Sodium thiosulphate (bleach neutralization)

Culture reagents

- Compressed carbon dioxide gas cylinder (pH control). The regulator supplied with each PBR will be a CGA320 threading in North America and BS341#8 threading internationally (see order confirmation). Ensure cylinder thread matches regulator supplied.
- Nutrient stock solutions. For marine algae, the commercially available (Proline or Fritz) Guillard's F/2 concentrates are suitable: Part A, B & silicate (diatoms only) are all required. If using other nutrient stock solutions, please contact support@industrialplankton.com with your recipe.
 - The PBR is equipped with two nutrient stock solution filters. The right nutrient pump is not filtered, as Silica will clog the 0.2 μm filter.
- Clean algae culture (inoculum)
 - Equip inoculum vessel with a Female Quick Disconnect Fitting (two supplied per PBR with ¾" ID barb) and 0.2 μm air inlet and vent filters.



Supplies and Equipment

	non-filtered nutrient stocks)
	ninum foil and/or autoclavable bags (maintaining sterility of autoclaved
	ponents)
	ay bottle(s) with adjustable mist setting (alcohol sanitizing fittings)
□ 5 ga	l/20 L bucket (calibrating Harvest Pump)
1 00	0 mL & 100 mL graduated cylinders (calibrating Nutrient Pumps)
☐ Soft	cotton rags (wiping acrylic)
□ 0 - 5	500 g capacity scale (weighing reagents)
☐ Carl	oon dioxide monitor/alarm (recommended around any CO ₂ source)
☐ Pall	et jack and/or forklift (moving/placing empty PBR)
☐ Step	pladder
Infract	ructure
☐ Stur	dy floor
	Floor must safely support at least working weight of PBR.
-	Sealed concrete is recommended for biosecurity.Floor should slope to drainage.
☐ Floo	
<u> </u>	 PBR is designed to gravity drain during cleaning.
	■ Drain should be within 7 m of the PBR.
	Cleaning fluids should drain to waste, not to waterways.
	Never allow cleaning fluids to mix.
	Most cleaners used are biodegradable through dilution, but check local
	discharge regulations to ensure compliance.
☐ Wire	eless Router with internet connection in range of PBR
	 Network firewalls may interfere with remote connection. Contact your
	network administrator to ensure remote access will be usable if desired.
☐ Prot	ection from the elements
•	PBRs should be protected from direct sunlight, or rain. Greenhouses are
	acceptable, but may require a larger integrated chiller or shadecloth to
□ \/a==	maintain acceptable culture temperatures.
	tilation PBRs should be well ventilated to disperse heat produced. Ventilation is
'	beneficial when working with compressed gases and cleaning agents.



- ☐ Continuous water supply for culture media (salty or fresh water)
 - Customer to equip supply with ½" female national pipe thread ball valve.
 - Water Supply should be constant (24 hours/day, 7 days/week)
 - Pressure: 0.3-4.1 bar (5-60 PSI)
 - Water chemistry needs to suitable for algae growth (free of chlorine or other chemicals which inhibit growth)
 - Water should be filtered down to 1µm (nominal) prior to header tank.
- ☐ Freshwater supply for cleaning
 - Free of particulate (<20 μm)
 - Municipal water is generally suitable
 - Equipped with a hose for filling PBR intermittently during cleaning process
- ☐ **Two (2)** Power Electrical Receptacles sized and installed in accordance with local regulations
 - In-use covers are recommended for all outlets installed around water sources. Ensure in-use covers are large enough to accommodate the plugs.
 - Main PBR Plug is equipped with GFCI circuit breaker inside control box. Please notify <u>support@industrialplankton.com</u> if facility breaker has GFCI installed and the internal GFCI will be removed.
 - 1 x Main PBR Power Receptacle (Refer to Order Confirmation)

Main PBR Options:	Nema L5-30	Nema L6-20	IEC 60309
Circuit Requirement:	110-120 V 60Hz Max Draw: 2690W / 115V / 23.6A	220V-240V (NOT 208V) 60Hz Max Draw: 2910W / 230V / 13.8A	220V 50Hz Max Draw: 2600W / 220V / 12.2A



■ 1 x Header Power Receptacle (Refer to Order Confirmation)

Header Outlet	Type B	Type F	Type G	Type I	Type M
Circuit	110-120 V	220-240 V	220-240 V	220-240 V	220-240 V
	60 hz	50 hz	50 hz	50 hz	50 hz
	Max Draw:	Max Draw:	Max Draw:	Max	Max Draw:
	450W /	450W / 230V	450W / 230V	Draw:450W /	450W / 230V
	115V / 4.3A)	/ 2.2A)	/ 2.2A	230V / 2.2A	/ 2.2A

■ Header does not include an internal GFCI. Ensure all wall receptacles are equipped with GFCI style outlets around liquids.

Specifications

Volume	115 L / 30 gal
Working Volume	100 L / 27 gal
Minimum Volume	16 L / 4.2 gal
Max Power Consumption	2,460 W
Max Heat Energy	8,400 BTU / hr
Avg. Steady State Power Req.	1,500 W
Avg. Heat Energy	5,100 BTU / hr
Available Voltages (AC)	110-120 / 220-240 V
Available Frequencies	60 / 50 Hz
Ideal Working Space (LxWxH)	2.03 m x 2.54 m x 1.98 m 80 in x 100 in x 78 in
Operational Weight	500 kg / 1,100 lb
Doorway Clearance (WxH) No dissasebly required	1.02 m x 1.79 m 40 in x 70 ¹ / ₂ in
Doorway Clearance (WxH) Some dissasebly required	0.87 m x 1.79 m 33 ⁷ / ₈ in x 70 ¹ / ₂ in
Doorway Clearance (WxH) Complex dissasebly required	0.76 m x 1.79 m 30 in x 70 ¹ / ₂ in
PBR 100L as shipped (LxWxH) (Wt)	1.70 m x 1.22 m x 2.08 m 67 in x 48 in x 82 in 363 kg / 800 lbs
Header Tank as shipped (LxWxH) (Wt)	1.19 m x 1.19 m x 1.80 m 47 in x 47 in x 71 in 140 kg / 300 lbs

Specifications subject to change without notice



Dimensions

